



The Imaging Phenomics  
Company®

# iBiopsy® AI-powered Lung Cancer Screening

Yan Liu, Chief Medical Officer

Benoit Huet, AI & Data Science Director, iBiopsy®

*Credits to the iBiopsy® AI/DS and Clinical Science teams*

## RSNA AI Theater



*“Many diseases, including cancer, do not kill if diagnosed early”*

We help conquer cancer and other health-threatening diseases through routine AI-powered digital biomarkers and imaging services for drug development

**Our Growth:** Powered by proprietary technology, strong KOL connections, and medical, scientific, technology partnerships

**Our People:** As of October 2021, 180+ employees worldwide (EU, US & China), 30% working in R&D

**Our locations:** HQ in France with subsidiaries in the US and in China

**Our 2 Business Units:**

- **iCRO:** image management for oncology trials
- **iBiopsy<sup>®</sup>:** AI-powered SaMD as novel disruptive digital biomarkers


**Our iBiopsy<sup>®</sup> clinical partners**

ASSISTANCE PUBLIQUE  HÔPITAUX DE PARIS

**UC San Diego**

# Solutions for Disease Diagnosis & Monitoring

## *Transforming the science of medical imaging*




# iBiopsy<sup>®</sup> Lung Cancer Screening

*A unique AI-powered end-to-end CAdE/CADx SaMD for early diagnosis*

## Intended use

iBiopsy<sup>®</sup> - Lung Cancer Screening is an AI powered Software as a Medical Device designed to


- Automatically **detect** and **characterize** pulmonary nodules from CT images
- For each detectable lesion, results including location of the lesion, lesion segmentation, and all nodule characteristics
- Provide a validated risk assessment of the probability of lung cancer for each detected nodule
- A per case score is also provided for the entire CT exam
- Guide a physician for patient management



# iBiopsy<sup>®</sup> Lung Cancer Screening

*A unique AI-powered end-to-end CAdE/CADx SaMD for early diagnosis*

## Processing Pipeline



## Baseline Model

- Nodule Regions Of Interest (ROI) are extracted through manual detection/segmentation by radiologists
- These ROI are then processed by a proprietary deep network trained in fully supervised fashion
- This deep network model (based on Convolutional Neural Nets and fully-connected layers) outputs a malignancy probability

# iBiopsy<sup>®</sup> Lung Cancer Screening

*Automatically detect, segment & characterize nodules*

## Lung Segmentation


- Deep Neural Network trained on LIDC/Luna16 dataset
- Independent verification on PleThora
- DICE "similarity coefficient"  $\geq .98$

## Nodule Detection/Segmentation

- Experiments ongoing using LIDC/Luna16 and NLST datasets

## Nodule Characterization

- Nodule-wise characterization on a subset of NLST



*Malignant*

*Benign*



# iBiopsy<sup>®</sup> Lung Cancer Screening: Cutting-edge performance with innovative deep learning algorithm


## Nodule Characterization (Nov 23<sup>rd</sup> Press Release)

### Nodule-wise characterization on NLST sub-set

- › Cohort of 1737 patients (16249 nodules)
- › Training set: 1239 patients (11676 nodules)
- › Test set: 498 patients [330 benign, 168 malignant] (4573 nodules)

### Deep Architecture with outstanding performance



- › **Test AUC = 0.99**
  - Youden Index Operating Point<sup>\*</sup>:  
Sensitivity = 95.3 %, Specificity = 96.2 %
  - Operating point **A**: Sensitivity 98.6% at Specificity 80.2%
  - Operating point **B**: Sensitivity 97.3% at Specificity 91.4%



# iBiopsy<sup>®</sup> Lung Cancer Screening

## Stage 1 (1A+1B) Lung Cancer characterization performance

Unprecedented performance beyond the state of the art



Test AUC = 0.984



- Youden Index Operating Point <sup>\*</sup>:  
Sensitivity = 89.2 %, Specificity = 96.8 %
- Operating point **A**: Sensitivity 98% at Specificity 80.2%
- Operating point **B**: Sensitivity 93.1% at Specificity 96.2%



# iBiopsy<sup>®</sup> Lung Cancer Screening

## Stage 1A Lung Cancer characterization performance

*Unprecedented performance beyond the state of the art*



**Test AUC = 0.982**

- Youden Index Operating Point<sup>\*</sup>:  
Sensitivity = 85.3 %, Specificity = 97.4 %
- Operating point **A**: Sensitivity 97.7% at Specificity 80.2%
- Operating point **B**: Sensitivity 92.1% at Specificity 96.2%



iBiopsy®



**iBiopsy®** Lung Cancer Screening is an AI-powered Software as a Medical Device that has demonstrated very promising results, that outperform the state of the art.



### **iBiopsy®** LCS benefits:

- Fully automatic end-to-end imaging workflow support,
- Early pulmonary nodules identification leading to early treatment,
- False negative and false positive reductions,
- Unnecessary invasive procedure & healthcare spending reduction.



**iBiopsy®** LCS end-to-end results are expected in the coming weeks



**iBiopsy®** currently focuses on Lung Cancer Screening (LCS), Liver cancer early diagnosis (HCC), Liver cancer (HCC) recurrence prediction and Non-Alcoholic Steato Hepatitis (NASH) diagnosis and scoring.



# Thank you!

## Meet us at booth #4849 AI Showcase | South Hall | Level 3

