

Median Technologies & iBiopsy[®]

March 2nd, 2023 ECR AI Theater

Yan Liu, Chief Medical Officer Benoit Huet, AI & Data Science Director, iBiopsy[®]

Early Diagnosis Saves Lives

We are helping conquer cancer and fibrotic diseases by extracting unvaluable clinical insights from patients' medical images.

median

Our strength is our people: 230+ highly qualified professionals in the US, Europe and China, 25+ nationalities.

Our growth is powered by proprietary AI, computer vision and signal processing technologies, strong KOL connections, and medical, scientific, technology partnerships.

With **iBiopsy**[®], our AI-powered imaging platform for the development of Software as Medical Device, we enable clinicians to diagnose patients earlier. We intend to launch our iBiopsy[®] Lung Cancer Screening SaMD in the US in 2024.

Our i**CRO** imaging solutions and advanced **Imaging Lab** offer help our 80+ biopharma clients drive their oncology clinical studies toward successful approval, using Ai-driven image insights.

iBiopsy[®] "shifting the diagnostic paradigm"

Lung Cancer Screening

I-ELCAP study showed a 92% survival rate at 15Y when diagnosed at stage 1 vs. 5% for stage 4 Lack of diagnosis accuracy is a major hurdle to screening adherence & programs implementation

Facts & Figures

- 1st cancer killer worldwide 18% of all 2020 cancer deaths
- 1.8M deaths in 2020
- 5-year overall survival rate 18%
- <25% cases in stage 1 when 5-y survival rate is 68%-92%</p>
- >40% cases in stage 4 when 5-y survival rate is <10%
- Rising frequency among never-smokers, 20% in the US & UK
- Only 870K screenings performed in the US in 2021 6% compliance



Why iBiopsy[®]?

Why is LDCT screening % so low in the high-risk populations ?



M median

Our Differentiators iBiopsy[®], our end-to-end CADe/CADx SaMD





iBiopsy[®] LCS Intended Use



AI/ML technology-based end-to-end CADe/CADx SaMD / MDSW





Intended use

iBiopsy LCS is an **AI/ML technology-based end-to-end CADe/CADx** Software as Medical Device intended to allow **early detection and characterization** of lung nodules and to **drive clinical management of patients**.

It is intended to aid the identification of tumor onset at earliest stage, allowing a better patient care while avoiding unnecessary tests and procedures and reducing healthcare costs.

Lung & Liver SaMD Workflow

Available as Saas or On-Premise – 4 steps workflow



M median

iBiopsy[®] LCS CADe/CADx Processing Methodology



State-of-the-art Artificial Intelligence & Data Science detection segmentation and characterization



iBiopsy[®] Lung Cancer Screening

Automatically segment left and right lungs

Lung Segmentation

- Dataset:
 - > Training and testing LIDC-IDRI/LUNA16
 - verified and reannotated manually
- Model:
 - > Ensemble of Deep Neural Networks
- Performance on test set (80 scans):> DICE = 0.991





iBiopsy[®] Lung Cancer Screening Lung nodule detection (CADe)

Nodule Detection

LIDC/IDRI Cohort (Luna16 Grand Challenge)

- Lung CT scans of 888 patients annotated by 4 experienced radiologists for a total of 1,186 nodules
- > Training set: 800 patients
- > Test set: 88 patients
- > training/testing using 10-fold cross-validation method

Encoder-Decoder Deep Architecture

- > Results: (obtained by aggregating the ten test folds)
- > Sensitivity of 94.9%

for a false positive rate of 1 per CT scan.





iBiopsy[®] Lung Cancer Screening

Lung nodule segmentation and measurements

Nodule Segmentation

- Dataset (LIDC/IDRI + NLST):
 - Training: 58594 nodules
 - Test: 13952 nodules
- Model: Ensemble of Deep Neural Networks

Nodule Measurements

- Volume
- Mean Diameter



iBiopsy[®] LCS CADe/CADx End to End Results



End-to-end Lung nodule detection & characterization with outstanding performance

Nodule detection and Characterization

- Dataset
 - Training (LIDC/IDRI + NLST):
 - > 7687 patients (225 malignant) 158686 Nodules
 - Test (NLST):
 - > 2176 Patients (independent from Train)
 - > 37257 Nodules (172 malignant)
- Model: Combination of Deep Neural Networks (2D + 3D), 3D-Morphological and Radiomics features (Patents Pending)
- Nodule Level Performance
 - AUC = 0.974
 - Max Youden Index Operating Point *:

Sensitivity = 96.5 % Specificity = 97.2 %



iBiopsy[®] LCS Morphological Feature

A key characterization feature for Lung Nodules

• Encoding local nodule curvature

• Effectively captures the shape irregularities (spiculations)

Conclusion & Future Work

iBiopsy[®] Lung Cancer Screening is an AI/ML tech-based SaMD that has shown very promising results, that outperform the state of the art.

iBiopsy[®] LCS benefits:

-End-to-end imaging workflow support,

- -Early pulmonary nodules identification leading to early treatment,
- -False positive and false negative reduction,
- -Unnecessary invasive procedure & healthcare spending reduction.

iBiopsy[®] currently focuses on Lung Cancer Screening (LCS) & Incidental pulmonary nodules (IPN), Liver cancer (HCC), and Non- Alcoholic Steato Hepatitis (NASH) fibrosis diagnosis and scoring.

iBiopsy[®] LCS 510(k) FDA clearance expected by first half 2024

iBiopsy

median

Our Core Values

Leading innovation with purpose

Combine the spirit of innovation with our passion and conviction to help cure cancer and other debilitating diseases.

Committing to quality in all we do

Be dedicated to quality in everything we do. Quality begins with us and we are committed to it.

Supporting our customers in achieving their goals Listen to the needs of our customers and help make their goals our goals through our innovation, imaging expertise, superior services, and quality solutions.

Putting the patient first

There is a person at the other end of the images we analyze who is counting on us to do everything we can to help make them healthier.

> ALMDT EURONEXT GROWTH