

# 274P - AI/ML-BASED LUNG CANCER DETECTION AND CHARACTERIZATION FOR LUNG CANCER SCREENING: RESULTS FROM THE REALITY STUDY ON EARLY-STAGE LUNG CANCER

Authors: R. OSAROGIAGBON<sup>1</sup>, A. VACHANI<sup>2</sup>, C. GOTERA<sup>3</sup>, L. SEIJO<sup>4</sup>, G. BASTARRIKA<sup>4</sup>, E. OSTRIN<sup>5</sup>, J. DENNISON<sup>5</sup>, C. VOYTON<sup>6</sup>, P. BAUDOT<sup>6</sup>, E. GEREMIA<sup>6</sup>, P. SIOT<sup>6</sup>, V. LE<sup>6</sup>, RENOUST<sup>6</sup>, B. HUET<sup>6</sup>, V. BOURDES<sup>6</sup>

[1] Baptist Cancer Center, Memphis, TN, [2] Penn Medicine, Philadelphia, PA, [3] Hospital Universitario Fundación Jiménez Díaz, Madrid, Spain, [4] Clínica Universidad de Navarra, Madrid, Spain, [5] MD Anderson Cancer Center, Houston, TX, [6] Median Technologies, Valbonne, France.

## BACKGROUND

AI-assisted radiological reading is entering clinical practice, potentially revolutionizing lesion characterization and early lung cancer detection. Tumor size and cancer stage at diagnosis are key determinants of survival: small, early-stage cancers are more responsive to treatment, with better prognoses. An AI/ML-Based computer assisted detection and diagnostic algorithm (CADE/CADx), was designed to assist the detection and characterization of suspicious lung nodules during lung cancer screening (LCS).

We present a subset of the pivotal REALITY Study results focusing on cancer stage and cell type.


## METHODS

Retrospective cohort study of individuals meeting USPSTF LCS criteria at five academic centers in Europe and the USA.


LDCT exams analyzed with the CADe/CADx; output compared to radiologist consensus. Reference standard proven via histopathology or  $\geq 12$  months stability.

AI/ML-Based CADe/CADx developed so that suspicious nodules presented on the results report should receive additional follow-up.

## FIGURES OF INTEREST



**Figure 1:** Presentation of the global workflow of the CADe/CADx for lesion prediction, with CADe 3D convolutional neural network (CNN) detection and lung segmentation on the left, and with nodule segmentation and CADx on the right. At the end a results report is generated.



**Figure 2:** Flowchart of participant selection and diagnostic performance of the CADe/CADx algorithm. Of 1,551 potentially eligible participants, 1,147 underwent the reference standard, with the device achieving an AUC-ROC of 0.904, sensitivity of 80.1%, and specificity of 86.6% at the maximum Youden Index (MYI).

## RESULTS

Category	Value
<b>Total Patients (Count)</b>	• 1147
<b>Patient Region (Count, % Total)</b>	• EU: 336, 29.29% • USA: 811, 70.71%
<b>Cancer Diagnosis (Count, % Total)</b>	• <u>Non-Cancer</u> : 805, 70.18% • <u>Cancer</u> : 342, 29.82%
<b>Average Age (Years, SD)</b>	• <u>Entire Cohort</u> : 64.94 ± 6.35 • <u>Non-Cancer</u> : 63.95 ± 6.30 • <u>Cancer</u> : 67.26 ± 5.83
<b>Nodule Size (mm, SD)</b>	• <u>Non-Cancer</u> : 6.37 ± 2.91 • <u>Cancer</u> : 14.91 ± 5.66
<b>Cancer Stage (Count)</b>	• <u>Stage I</u> : 219 • <u>Stage II-IV</u> : 56 • <u>Not Available</u> : 67
<b>Cancer T-category (Count)</b>	• <u>T1</u> : 227 • <u>T2-4</u> : 43 • <u>Not Available</u> : 72
<b>Lung Cancer Cell Type (Count)</b>	• <u>NSCLC</u> : 213 • <u>Other Cell Type (SCLC Included)</u> : 64 • <u>Not Available</u> : 65

**Table 1:** Study Participants.

Lung Cancer Cell Type	AUC	Sensitivity (Cancer Recall Rate)
<b>NSCLC</b>	0.908	97.2%
<b>Other Lung Cancer Cell Types</b>	0.853	89.1%

**Table 2:** AUC for correct lesion characterization and cancer recall rate (sensitivity) according to cancer cell type.

Class	Sub-Class	Sensitivity (Cancer Recall Rate)	False Positive Detections Per Exam
<b>Cancer Stage</b>	Stage I	96.8%	0.80 FP/ Scan
	Stage II-IV	98.2%	0.86 FP/Scan
<b>Cancer T-Category</b>	T1	96.5%	0.80 FP/Scan
	T2-4	100%	0.86 FP/Scan

**Table 3:** Cancer recall rate and false positive marks per examination for patients according to cancer stage and cancer t-category

## CONCLUSIONS

The AI/ML-Based CADe/CADx showed high performance in detecting early lung cancer, high sensitivity for cancers of different stage, T-category and cell type, excellent cancer recall rate (>96.5%) on all cancer stages and t-categories. This CADe/CADx can assist in the accurate follow-up of suspicious lung nodules, optimizing the management of patients in lung cancer screening.

## CONFLICTS OF INTEREST

Presenter R. OSAROGIAGBON: consultancies and advisory Median Technologies.